

## MAMMOMAT Novation

**SP**

### Installation Protocol

System

Service provider Installation Protocol

© Siemens AG 2004

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights, including rights created by patent grant or registration of a utility model or design, are reserved.

**Document revision level**

The document corresponds to the version/revision level effective at the time of system delivery. Revisions to hardcopy documentation are not automatically distributed.

Please contact your local Siemens office to order current revision levels.

**Disclaimer**

The installation and service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such services.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are directed to contact one of the local offices of Siemens or one of its affiliates before attempting installation or service procedures.

1	<b>Scope of services</b>	4
	Delivery, installation, cabling and start-up	4
	Room preparation	4
	Delivery	4
	Installation	4
	System cabling	5
	Other work	5
2	<b>Installation Protocol</b>	6
	Site data	6
	Installation steps performed	7
	Activity report	9
	General remarks	10
3	<b>Changes to previous version</b>	11

## Delivery, installation, cabling and start-up

### Room preparation

- The following topics must be clarified with the project manager before the delivery date:

Building inspection by the service provider

Floor loading

Floor condition

Barriers, crane, tools.

- The project manager must be on site for information on arrival of the service provider.
- The orientation point is marked by the project manager.
- The current project plan is handed over to the service provider on site.
- Line voltage is available at the on-site system distributor.
- The system emergency off is installed by the customer.
- Cable ducts etc. must be installed.
- All room preparations by the customer are completed and the room is swept clean.
- Radiation protection partition is available.
- To be able to install the covers, the system must be started up partially. This requires a Siemens service engineer who connects and switches on the voltage to the MAMMO-MAT Novation<sup>DR</sup>.
- Seal the joints of the unit base with silicone.

<b>NOTE</b>
-------------

---

**All work must be performed in accordance with the supplied documents.**

---

### Delivery

- A long elevating truck is required for transporting the unit carrier on the pallet.
- The system is brought into the room.
- Unpacking and bringing in all system components of the unit.
- Inspecting the components for damage.

### Installation

- Taking along all required tools, measuring and installation aids according to the tool list.
- Positioning and aligning all system components according to the project plan or the statements of the PM.
- Installing accessories and options.

- Aligning, drilling and anchoring the unit base, if necessary, screw tight and tighten with the torque stated in the MAMMOMAT Novation<sup>DR</sup> installation instructions.
- Remove all transport protection devices.
- All fastenings must be tightened with the torques stated in the installation instructions and secured where necessary with Locite 221.
- All values must be entered according to instructions in the corresponding certificates.
- Install the detector

## System cabling

- Perform the cabling of the system components of the MAMMOMAT Novation<sup>DR</sup>, options, control consoles and monitors corresponding to the installation instructions.
- Connect foot switches.

## Other work

- Close all cable ducts.
- Clearing and cleaning the installation room (swept clean).
- Unpack accessories, if the premises allow this.
- Return all installation racks and packaging materials.
- Damage caused to the building by the work must be reported to the project manager/person responsible for the project.

**Site data**

SAP No.: \_\_\_\_\_ System Ser. No.: \_\_\_\_\_  
Customer: \_\_\_\_\_  
Country/city: \_\_\_\_\_  
Installation team: \_\_\_\_\_  
Project manager: \_\_\_\_\_  
Installation start: \_\_\_\_\_ Installation end: \_\_\_\_\_

Installation steps performed

	Work performed	
	yes	no
The system was installed according to the supplied documents.	<input type="checkbox"/>	<input type="checkbox"/>
Delivery of the system:		
All components were delivered. ....	<input type="checkbox"/>	<input type="checkbox"/>
The system delivery is complete, according to dispatch note. ....	<input type="checkbox"/>	<input type="checkbox"/>
All parts are undamaged. ....	<input type="checkbox"/>	<input type="checkbox"/>
The system emergency off is installed by the customer. ....	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>	<input type="checkbox"/>

	Work performed	
	yes	no
<b>Installation:</b>		
The stand is anchored to the floor. ....	<input type="checkbox"/>	<input type="checkbox"/>
The radiation shield is anchored to the floor. ....	<input type="checkbox"/>	<input type="checkbox"/>
Missing parts are noted in the installation protocol. ....		
The project plan was available. ....	<input type="checkbox"/>	<input type="checkbox"/>
The radiation shield and AWS are installed according to project plan . . . . .	<input type="checkbox"/>	<input type="checkbox"/>
Covers are mounted. ....	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>	<input type="checkbox"/>
<b>Cabling:</b>		
System cables are run and connected according to the documents. ....	<input type="checkbox"/>	<input type="checkbox"/>
All options are positioned and cabled . . . . .	<input type="checkbox"/>	<input type="checkbox"/>
Enough cable ducts were available and with a sufficient length . . . . .	<input type="checkbox"/>	<input type="checkbox"/>



## Activity report

Day	Start	End	E	M	V	R	KA	Work performed	Employees/ number
<b>Total (Mh)</b>									

### Explanation of the activity report

<b>Day / Start / End</b>	Date / start of work / end of work (time)
<b>E</b>	Bringing in the entire system in hours per/day
<b>M</b>	Installation of the entire system in hours per/day
<b>V</b>	Cabling of the system in hours per/day
<b>R</b>	Travelling time in hours per/day
<b>KA</b>	No order-related working time in hours per/day
<b>Work performed</b>	Work performed
<b>Employees/number</b>	Number of employees who were employed in the period
<b>Total (Mh)</b>	Man hours, number of accumulated hours multiplied by number of employees

General remarks

A copy of this remains at the system and is filed in the Assembly/Installation folder.

Place, date	Installation team employee	Siemens project manager/ person responsible for the project
	Name in block letters	Name in block letters

Chapter	Page	Change
All	All	Converted to DMS

